



**VESALIO**<sup>TM</sup>

*is dedicated to*

**CHANGING OUTCOMES, CHANGING LIVES**

*by advancing the care of patients suffering from vascular occlusion by providing physicians superior technology designed to improve clinical outcomes*

**enviast**

**THE ART OF CAPTURING ALL CLOT TYPES<sup>2</sup>**

# THE NEED



Large thrombus burden (LTB) is a **common** and **critical** occurrence in Acute Coronary Syndrome (ACS)

**60%**

of patients with ACS<sup>1</sup>

**75%**

of patients presenting  
with STEMI  
(ST-elevation myocardial  
infarction)

**x2**

higher risk of  
mortality

**x2 - x4**

higher risk of  
major adverse  
events<sup>2</sup>

**Studies show thrombus aspiration alone may increase stroke re-occurrence<sup>3</sup>**

**There is no established technique for managing LTB in ACS**

1. Napodano M, Dariol G, Al Mamary AH, Marra MP, Tarantini G, D'Amico G, Frigo AC, Buja P, Razzolini R, Iliceto S. Thrombus burden and myocardial damage during primary percutaneous coronary intervention. Am J Cardiol 2014;**113**(9):1449-56.
2. Singh M, Berger PB, Ting HH, Rihal CS, Wilson SH, Lennon RJ, Reeder GS, Bresnahan JF, Holmes DR, Jr. Influence of coronary thrombus on outcome of percutaneous coronary angioplasty in the current era (the Mayo Clinic experience). Am J Cardiol 2001;**88**(10):1091-6
3. Neumann FJ, Sousa-Uva M, Ahlsson A, Alfonso F, Banning AP, Benedetto U, Byrne RA, Collet JP, Falk V, Head SJ, Juni P, Kastrati A, Koller A, Kristensen SD, Niebauer J, Richter DJ, Seferovic PM, Sibbing D, Stefanini GG, Windecker S, Yadav R, Zembala MO, Group ESCSD. 2018 ESC/EACTS Guidelines on myocardial revascularization. Eur Heart J 2019;**40**(2):87- 165.

## 1. TREAT ALL LTB OCCLUSIONS

FROM SOFT CLOTS  
THAT EASILY DISINTEGRATE  
TO HARD, FIBRIN-RICH  
CLOTS THAT CANNOT BE  
INGESTED

## 2. IMPROVE PROCEDURAL PERFORMANCE

HIGHER 1ST PASS SUCCESS  
BETTER TIME TO  
RECANALIZATION

## 3. PROVIDE EASE OF USE

REAL TIME FEEDBACK  
DURING RETRIEVAL  
SYNERGISTIC WITH  
ASPIRATION

**To achieve better patient outcomes**



en<sup>v</sup>ast

**DROP ZONE™**  
**THE CLOT INSIDE**



**CE approved** for temporary endovascular use to restore blood flow in patients experiencing thrombosis symptoms in the coronary vasculature (December 2019)



**Uniquely designed  
to capture ALL CLOT TYPES  
INSIDE the device structure**


CONFIDENTIAL, FOR PRESENTATION PURPOSES ONLY. DO NOT DISTRIBUTE.

### FLOW RESTORATION ZONE

Clot initiation zone with optimized radial force & large cell size

### DROP ZONES

Entry points for organized/ hard thrombi



### BALANCED DESIGN

Optimized radial force balanced with large openings & closed ends

### SMART MARKERS

2 per Drop Zone, for real-time feedback during retrieval

### CLOSED DISTAL TIP

Clot gets inside, clot stays inside!

## CE approved

The Vesalio enVast Mechanical Thrombectomy System is indicated for endovascular temporary use to restore blood flow in patients who are experiencing symptoms of thrombosis in the coronary vasculature.



# enVast PORTFOLIO




Product Number	Product Name	Labeled Device Diameter (mm)	Labeled Device Length (mm)	Self Expanded Device Diameter (mm)	Recommended Vessel Diameter (mm)	Pusher Length	Introducer Microcatheter Minimum Inner Diameter
EV-4537-F2RR	enVast 4.5 x 37 mm	4.5	37	4.5	$\geq 2.0$ and $\leq 4.5$	180cm	.021"
EV-4546-F3RR	enVast 4.5 x 46 mm	4.5	46	4.5	$\geq 2.0$ and $\leq 4.5$	180cm	.021"
EV-4030-F2RR	enVast 4.0 x 30 mm	4.0	30	4.0	$\geq 2.0$ and $\leq 3.5$	180cm	.021"
EV-4038-F3RR	enVast 4.0 x 38 mm	4.0	38	4.0	$\geq 2.0$ and $\leq 3.5$	180cm	.021"
EV-6035-F2RR	enVast 6.0 x 35 mm	6.0	35	6.0	$\geq 3.5$ and $\leq 6.0$	180cm	.027"
EV-6044-F3RR	enVast 6.0 x 44 mm	6.0	44	6.0	$\geq 3.5$ and $\leq 6.0$	180cm	.027"

# enVast PORTFOLIO



enVast 4.0 x 30 mm, 2 Drop Zones, Full length: 39 mm




Vessel diameters 2.0 – 3.0 mm

Ideal for side branches or distal tracts of main vessels

Compatible with: 0.021" MC

enVast 4.5 x 37 mm, 2 Drop Zones, Full length: 57 mm



Vessel diameters 2.0 – 4.5 mm

Ideal for proximal to mid tracts of main branches

Compatible with: 0.021" MC

enVast 4.5 x 46 mm, 3 Drop Zones, Full length: 66 mm




Vessel diameters 2.0 – 4.5 mm

Ideal for proximal to mid tracts of main branches  
if thrombus burden is very large or difficult to quantify

Compatible with: 0.021" MC

enVast 6.0 x 35 mm, 2 Drop Zones, Full length: 55 mm



Vessel diameters 3.5 – 6.0 mm

Ideal for large RCA, left main or ectatic vessels  
or every time the vessel RVs is greater than 4.00

Compatible with: 0.027" MC



# CURRENT STATUS



## First in man (FIM) completed

- 61-patient case series submitted for publication

## FIM presented

- at the Transcatheter Cardiovascular Therapeutics (TCT) Meeting in November

Next Step: Starting European Registry



## FIRST-IN-MAN

- Two tertiary centers in Switzerland (Bern, Lugano)
- 61 consecutive ACS patients with LTB (TTG  $\geq$  3)

## EFFICACY ENDPOINTS

- ST-segment elevation resolution and TIMI flow
- TIMI Thrombus Grade and Myocardial Blush Grade by an independent core laboratory (MedsStar Washington Hospital Centre, Washington DC, USA)




## SAFETY ENDPOINTS

- device and procedure-related adverse events
- MACCE and bleedings at 30 days

# FIRST-IN-MAN EFFICACY OUTCOMES




enVast deployment → immediate reperfusion in 85% of cases



# FIRST-IN-MAN EFFICACY OUTCOMES



**ST resolution > 50% in 71.7% of cases**  
**> 70% in 43.5% of cases**



# FIRST-IN-MAN SAFETY OUTCOMES



- Cardiovascular death in two patients (3.3%) - in cardiogenic shock at admission;
- No major procedure-related adverse events (coronary dissection, coronary perforation, cardiac tamponade, coronary occlusion, life threatening arrhythmias)
- 14 (23%) non-flow-limiting coronary spasms, resolved with intracoronary nitrates
- 1 (1.6%) unplanned revascularization at 30 days (stent under-expanded);
- 1 (1.6%) case (without continuous aspiration) of side-branch embolization requiring additional stent retrieval with complete vessel reperfusion;
- 1 (1.6%) transient ischemic attack at day 29, after a conventional staged PCI

<b>Procedural outcomes</b>	<b>n=61</b>
Coronary dissection	0
Coronary perforation	0
Coronary occlusion	0
Coronary spasm	14 (23.0)
Flow limiting spasm	n=14, 0 (0)
Spasm resolution	n=14, 14 (100)
Embolization	1 (1.6)
Embolization resolution	n=1, 1 (100)
Cardiac tamponade	0
Life-threatening arrhythmias needing treatment	0

<b>Clinical Outcomes</b>	<b>n=61</b>
Death	2 (3.3)
Cardiovascular death	2 (3.3)
Non-cardiovascular death	0
Myocardial infarction	0
Unplanned revascularisation (any)	1 (1.6)
Definite Stent thrombosis	0
Cerebrovascular events	1 (1.6)
Stroke (any)	0 (0)
Transient ischemic attack	1 (1.6)
Bleeding BARC 3 or 5	0
Bleeding BARC 2	3 (4.9)
access site	3 (4.9)
non-access site	0

Data are presented as absolute numbers (percentage).

# REGISTRY CONSIDERATIONS

en<sup>v</sup>ast THROMBECTOMY DEVICE REGISTRY STUDY IN ST-SEGMENT ELEVATION MYOCARDIAL INFARCTION



A retrospective, multi-center registry study designed to assess the safety and effectiveness of the en<sup>v</sup>ast coronary thrombectomy system as an adjunctive measure to conventional intervention in subjects presenting with ST-segment elevation myocardial infarction (STEMI)

Up to 200 subjects at up to 15 sites

en<sup>v</sup>ast will be deployed as the first measure to obtain reperfusion at the occlusion site up to 3 times

± Conventional Treatment  
(ballooning, manual aspiration thrombectomy, stenting)

# REGISTRY CONSIDERATIONS

enVast THROMBECTOMY DEVICE REGISTRY STUDY IN ST-SEGMENT ELEVATION MYOCARDIAL INFARCTION



## Identification of Candidate Sites

- Hospital name & address
- Physician Name & contact information (principal investigator)
- Number of interventionalists that would contribute to the study at site
- STEMI case load/ year
- How do they currently treat large clot burden STEMI cases?
- Why should we consider this site?



## Clinical Team to Present Protocol to the Site PI & Team

- Can they work with the protocol?
- Can they get informed patient consent within 72hrs from case?
- Documents needed for approval?
- Estimated timeline to approval?
- Buying enVast: any hurdles to consider?



## Training / Evaluation

- Expect 3-5 cases/ physician to overcome the learning curve
- Cases can be done while the site is awaiting approval of the protocol

# DISCUSSION OF NEXT STEPS



Selling – what needs to be considered?

Identification of sites: please come back with list (excel file will be provided after the call)



# NATURE STUDY SITES TO EXCLUDE FROM OPTIONS FOR THE REGISTRY



Country	City	Hospital	PI
Switzerland	Bern	University Hospital Inselspital of Bern	Jonas Häener
	Lugano	Cardiocentro Ticino	Marco Valgamigili
	Genève	Hôpitaux Universitaires de Genève (HUG)	Juan Fernando Iglesias
	Luzern	Lucerne Cantonal Hospital	Florim Cuculi
Italy	Genova	IRCCS Ospedale Policlinico San Martino (referred to as "Policlinico")	Crimi
	Milan	Niguarda Hospital	Jacopo Oreglia
	Milan	Humanitas Hospital	Antonio Colombo
	CONA (FE)	Azienda Ospedaliero Universitaria di Ferrara	Matteo Tebaldi

