

SCHEDA TECNICA FONTE DI LUCE LED
Modello DAYMARK

| | |
|------------------------------|--|
| Descrizione: | Lampada frontale a LED portatile ad alta intensità |
| Ditta Produttrice: | BFW INC. U.S.A. |
| Ditta Importatrice: | Medical Instruments S.p.A. |
| Modello: | Daymark – Lampada frontale a LED portatile ad alta intensità |
| Codice: | BFW 9070 (Sistema Daymark con Caschetto Fiji completa di n. 2 batterie e carica batterie) |
| Classe di protezione e tipo: | Classe I |
| Applicazioni Direttive EEC: | 93/42/EEC |
| Codice CND: | Z129004 |
| Registrazione Ministero RDM: | 2062055 |

CARATTERISTICHE TECNICHE:

| | |
|-------------------------------|--|
| Vita lampada: | 30.000 ore |
| Temperatura di colore: | 5.000° K |
| Intensità luminosa: | da 60.000 lux a 40 cm a 260.000 lux a 25 cm |
| Diametro Spot: | diametro regolabile da 2.5 a 10.16 cm a una distanza di 40 cm |
| Controllo Intensità Luminosa: | 5 livelli di intensità impostati |
| Batteria: | Ioni di litio, gruppo a 4 celle |
| Durata Batteria: | 5 ore alla massima intensità |
| Tempo di carica Batteria: | 2-3 ore per la carica completa |
| Cicli di ricarica: | più di 300 |


DIMENSIONI E PESO:

Peso: **Batteria 226 gr**
Ottica 141 gr

Garanzia: **12 mesi**

DESTINAZIONE D'USO:

La Fonte di luce a LED è indicata per fornire illuminazione per la visione durante interventi chirurgici, diagnostici o terapeutici o ambulatoriali.



DAYMARK™ | SURGICAL

PORTABLE LED HEADLIGHT SYSTEM



SEE BETTER
SEE FARTHER
SEE WHY



1971
2021

CELEBRATING
FIFTY YEARS OF
ILLUMINATION


SINCE 1971
ENHANCING
YOUR VISION

Powered by
Long Island
Technology
Group

DAYMARK™ SURGICAL

HIGH-INTENSITY, PORTABLE LED SURGICAL HEADLIGHT WITH ADJUSTABLE SPOT, POWERED BY A 4-CELL BATTERY

- Illumination Range: ~60,000 LUX at 16" (40 cm) to ~260,000 LUX at 10" (25 cm)
- Adjustable Spot Diameter: 1"-4" (2.5 - 10.16 cm) at a working distance of 16" (40 cm)
- Color Rendering Index: 78 CRI
- Color Temperature: 5,000 K
- Intensity Settings: 5 Levels
- 4-Cell Battery
- Battery Duration at Full Intensity: 5 hours
- Battery Expected Life Span: > 300 Charges
- LED Expected Life Span: > 10,000 hours
- Optic Weight: 5 oz (141 g)
- Battery Weight: 8 oz (226 g)



DAYMARK PORTABLE LED HEADLIGHT SYSTEM INCLUDES:

- Daymark Surgical Headlight on Fiji Headband
- Two 4-Cell Batteries
- Charging Cradle, Power Supply
- Holster, Linkage Kit

DAYMARK PORTABLE LED HEADLIGHT SYSTEM PRODUCT CODES:

- | | |
|-------------|---|
| BFW 9070 | Daymark Headlight System with Fiji Headband with Two 4-Cell Batteries |
| BFW 5440 | Daymark Protective Carrying Case, molded interior provides space for all components |
| BFW 9022ACC | Daymark Mini Battery Charger |
| BFW 9015 | Daymark 4-Cell Battery |
| BFW 9076 | Daymark Holster |

EXCELLENT FOR USE IN THE FOLLOWING DEPARTMENTS

HOSPITAL

- ENT
- GENERAL SURGERY
- CARDIOSURGERY
- NEUROSURGERY
- ORTHOPEDIC SURGERY
- PEDIATRIC SURGERY
- SPINE SURGERY
- THORACIC SURGERY
- TRANSPLANT SURGERY
- VASCULAR SURGERY

SURGERY CENTER

- ORTHOPEDIC SURGERY
- PLASTIC SURGERY

OFFICE

- DENTAL EXAM
- ORAL/MAXILLOFACIAL



BFWinc.com

800-717-4673 | +1-502-899-1808

445 BAXTER AVENUE | SUITE 175

LOUISVILLE, KY 40204 USA



DAYMARK TESTED AND CERTIFIED TO:
IEC 60601-1 Edition 3, IEC 60601-1-2 Edition 4
ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10),
CAN/CSA-C22.2 NO 60601-1- (2008), CE

ENGINEERED AND ASSEMBLED IN THE U.S.A. © 2020 BFW, Inc.
ALL TECHNICAL SPECIFICATIONS ARE SUBJECT TO CHANGE WITHOUT NOTICE