

"A REAL ADDED VALUE TO
YOUR INTERVENTIONAL
ARMENTARIUM."



Pallone per PTCA per lesioni complesse

- Basso profilo, rivestimento idrofilico del pallone per lesioni complesse
- RBP 21 atm
- 1 singolo marker centrale

- a) Profilo d'ingresso: 0.016"
- b) Crossing profile: 0.0205"
- c) Rivestimento idrofilico sul pallone e sullo shaft distale

SIS  MEDICAL
Swiss | Interventional | Systems

NIC 1.1[®] *hydro*

Ordering Information
NIC 1.1[®] hydro CTO-Balloons (RX)

Article Number	Balloon Diameter	Length	RBP	Marker	Min. Guiding
110-006-134	11 mm	6 mm	21 atm	1 single middle	5F
110-010-134	11 mm	10 mm	21 atm	1 single middle	5F
110-015-134	11 mm	15 mm	21 atm	1 single middle	5F
110-020-134	11 mm	20 mm	21 atm	1 single middle	5F



SIS MEDICAL

Swiss | Interventional | Systems

Hungerbühlstrasse 12a | CH-8500 Frauenfeld
Tel. +41 (0)52 245 09 90 | Fax +41 (0)52 245 09 91
info@sis-medical.com | www.sis-medical.com

Copyright © by SIS Medical AG
Currently not available for sale in USA.
Subject to alteration due to technical advances.
All rights reserved.
All products are developed and produced by
SIS Medical AG.